FoundationOne CDx版 C-CAT調查結果version1.0

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2019/7/23

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The screenshot displays a FoundationOne CDx report. At the top, it shows the patient's name (Chugai Unique ID), tumor type (Lung adenocarcinoma), and report date (01 Mar 2019). Below this, there are sections for Patient, Physician, and Specimen information, all of which are marked as 'Not Given'. The main section is titled 'Companion Diagnostic (CDx) Associated Findings' and is divided into two columns: 'GENOMIC FINDINGS DETECTED' and 'APPROVED THERAPEUTIC OPTIONS IN JAPAN'. Under 'GENOMIC FINDINGS DETECTED', there is one entry: 'EGFR L858R'. Under 'APPROVED THERAPEUTIC OPTIONS IN JAPAN', there are two entries: 'Afatinib maleate' and 'Erlotinib hydrochloride'. Below this, there is a section titled 'OTHER ALTERATIONS & BIOMARKERS IDENTIFIED' which contains a list of findings: 'Microsatellite Status MS-Stable', 'Tumor Mutational Burden 11 Muts/Mb', 'CDKN2A/B loss', 'EGFR amplification', 'PTCH1 T416S', 'RBM10 Q494*', and 'TP53 R267P'. A disclaimer at the bottom states: 'Refer to appendix for limitation statements related to detection of any copy number alterations, gene rearrangements, MSI or TMB result in this section. Please refer to appendix for Explanation of Clinical Significance Classification and for variants of unknown significance (VUS).'

GENOMIC FINDINGS DETECTED	APPROVED THERAPEUTIC OPTIONS IN JAPAN
EGFR L858R	Afatinib maleate Erlotinib hydrochloride

OTHER ALTERATIONS & BIOMARKERS IDENTIFIED
Results reported in this section are not prescriptive or conclusive for labeled use of any specific therapeutic product. See professional services section for additional information.

Microsatellite Status MS-Stable [§]	PTCH1 T416S
Tumor Mutational Burden 11 Muts/Mb [§]	RBM10 Q494*
CDKN2A/B loss [§]	TP53 R267P
EGFR amplification [§]	

[§] Refer to appendix for limitation statements related to detection of any copy number alterations, gene rearrangements, MSI or TMB result in this section.
* Please refer to appendix for Explanation of Clinical Significance Classification and for variants of unknown significance (VUS).

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