# PHILIPS

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Engineer Name: YOJIRO MOTOYAMA

Device Status: Pass

Event Type: Predictive / Preventative Maintenance

Service Activity: Preventative Maintenance

### EQUIPMENT DETAILS

Product Name: DigitalDiagnost 4 High Performance

Product Number: 712031

Serial Number: SN16000415

Installed Product: 70301047

UDI:

	Min	Max	UoM	Measurement	Result	Comments
PM Start						
PM Start						
Functional Status					PASS	
PM Table TH						
Table TH (24 Month)						
Functional Status					PASS	
Table Base					PASS	
Bucky Unit					PASS	
Final Steps					PASS	
PM Wall Stand VS						
Wall Stand VS (12 month)						
Fixations					PASS	
Wall Stand VS (24 month)						
Functional Status					PASS	
Covers					PASS	
Guide Rails and Castors					PASS	
Vertical Collision Check					PASS	
Bucky Unit					PASS	
Wrap Up					PASS	

	Min	Max	UoM	Measurement	Result	Comments
<b>PM Ceiling Suspension</b>						
CSM/CSM3 (12 month)						
Functional status					PASS	
Covers					PASS	
Fixations					PASS	
Wrap up					PASS	
CSM/CSM3 - Castors (Visits 1, 3, 5, 7, 9)						
Castors					PASS	
<b>PM Workstation</b>						
AWS - Acquisition Work Spot						
Functional Status					PASS	
Philips Remote Service Functionality Checked					PASS	
Automatic Export of Log File Configured					Yes	
Cleaning AWS and Accessories					PASS	
<b>PM M-Cabinet</b>						
M-Cabinet						
Functional Status					PASS	
Cable connections					PASS	
Tube calibration					PASS	
Backup of Generator data					PASS	
Final steps					PASS	
<b>PM X-ray Tubes</b>						
SRO/RO/SRM X-ray Tube Assemblies						
Preparatory Work					PASS	
Mechanical Status X-Ray Tube Assy.					PASS	
Status of Cables and Connectors					PASS	
Final Work					PASS	
<b>PM Collimators NICOL/Galileo/XD</b>						
Collimators NICOL/Galileo/XD (12 month)						
Collimator Lamp Replacement					PASS	
Safety and Radiation Safety Check					PASS	
Functional Check					PASS	

	Min	Max	UoM	Measurement	Result	Comments
PM Collimators NICOL/Galileo/XD						
Collimators NICOL/Galileo/XD (24 month)						
General Inspection					PASS	
Cleaning and Lubrication					PASS	
IQA Table TH						
IQA Checks Table						
X-ray Geometry Test					PASS	
Detector Performance Fixed Detector					PASS	
Amplimat Dose Test Table						
Amplimat dose test - Grid					110/8/40	
Amplimat dose test - Attenuation factor				1.9		
Amplimat dose test - S200			μGy	9.61		
Amplimat dose test - S400			μGy	4.57		
Amplimat dose test - S800			μGy	2.45		
Amplimat dose test					PASS	
IQA Wall Stand						
IQA Checks Wall Stand						
X-ray Geometry Test					PASS	
Detector Performance Fixed Detector					PASS	
Amplimat Dose Test Wall Stand						
Amplimat dose test - Grid					180/12/40	
Amplimat dose test - Attenuation factor				1.7		
Amplimat dose test - S200			μGy	8.34		
Amplimat dose test - S400			μGy	4.07		
Amplimat dose test - S800			μGy	1.95		
Amplimat dose test					PASS	
IQA Other Components						
IQA Checks Other Components						
Monitor Performance					PASS	
DAP Indication					PASS	
Check accuracy of loading factors					PASS	

	Min	Max	UoM	Measurement	Result	Comments
ES Table TH						
Protective Earth Table TH						
2) Potential Balance Point SX80		300	mOhm	55	Pass	
4) Foot Switch		300	mOhm	67	Pass	
5) Transversal carriage at SAYX 13		300	mOhm	66	Pass	
7) Service Cover		300	mOhm	21	Pass	
Protective Earth Bucky unit Table (Fixed Detector)						
2) Grid ejector		300	mOhm	22	Pass	
5) Rear cover screw		300	mOhm	32	Pass	
Touch Current Table TH						
2) Potential Balance Point SX80		20	μA	0.1	Pass	
3) Table Top Side Rail (Front)		20	μA	0.1	Pass	
3) Table Top Side Rail (Back)		20	μA	0.1	Pass	
4) Foot Switch		20	μA	0.1	Pass	
5) Transversal carriage at SAYX 13		20	μA	0.1	Pass	
7) Service Cover		20	μA	0.2	Pass	
Touch Current Bucky unit Table (Fixed Detector)						
2) Grid ejector		20	μA	0.2	Pass	
5) Rear cover screw		20	μA	0.1	Pass	
ES Wall Stand VS						
Protective Earth Wall Stand VS						
2) Column screw		300	mOhm	55	Pass	
3) Wall box cover screw		300	mOhm	45	Pass	
Touch Current Wallstand VS						
2) Column screw		20	μA	0.1	Pass	
3) Wall box cover screw		20	μA	0.1	Pass	
4) Rear cover screw wall stand display		20	μA	0.1	Pass	
5) Holder screw remote control station		20	μA	0.1	Pass	
ES Bucky Units Wall Stands						
Protective Earth Bucky unit Wall Stand (Fixed Detector)						
2) Grid ejector		300	mOhm	33	Pass	
6) Rear cover screw		300	mOhm	54	Pass	

	Min	Max	UoM	Measurement	Result	Comments
<b>ES Bucky Units Wall Stands</b>						
Touch Current Bucky unit Wall Stand (Fixed Detector)						
2) Grid ejector		20	μA	0.2	Pass	
6) Rear cover screw		20	μA	0.2	Pass	
<b>ES Ceiling Suspension</b>						
Protective Earth CS						
4) Union nut of the high voltage cable		300	mOhm	33	Pass	
6) Accessory rail screw		300	mOhm	33	Pass	
7) Rear screw of the tube control handle		300	mOhm	41	Pass	
Touch Current CS						
4) Union nut of the high voltage cable		20	μA	0.1	Pass	
5) Key switch		20	μA	0.1	Pass	
6) Accessory rail screw		20	μA	0.1	Pass	
7) Rear screw of the tube control handle		20	μA	0.1	Pass	
8) Clutch disk		20	μA	0.2	Pass	
<b>ES M-Cabinet</b>						
Protective Earth M-Cabinet CXA						
2/B) Front cover		300	mOhm	66	Pass	
3/C) Top cover		300	mOhm	67	Pass	
4/D) Rear cover		300	mOhm	76	Pass	
5/E) Left cover		300	mOhm	76	Pass	
6/F) Right cover		300	mOhm	60	Pass	
Touch Current M-Cabinet CXA						
2/B) Front cover		20	μA	0.1	Pass	
3/C) Top cover		20	μA	0.1	Pass	
4/D) Rear cover		20	μA	0.2	Pass	
5/E) Left cover		20	μA	0.2	Pass	
6/F) Right cover		20	μA	0.2	Pass	
<b>ES Workstation</b>						
Protective Earth Acquisition Workspot						
2) Screw (rear side of PC)		300	mOhm	55	Pass	

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

## ES Wall Connection Box

### Protective Earth Wall Connection Box

2) Upper Left Bezel		300	mOhm	61	Pass	
3) Lower Left Bezel		300	mOhm	62	Pass	
4) Upper Right Bezel		300	mOhm	56	Pass	
5) Lower Right Bezel		300	mOhm	66	Pass	
6) Front Cover Screw		300	mOhm	55	Pass	

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