PHILIPS

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Engineer Name: TOMOKI ABE

Device Status: Pass

Event Type: Predictive / Preventative Maintenance

Service Activity: Visit 7a

PM schedule for this site: PM Plan Template DigitalDiagnost 4.x- High- 20 visits / 120 months

EQUIPMENT DETAILS

Product Name: DigitalDiagnost 4 High Performance

Product Number: 712031

Serial Number: SN16000415

Installed Product: 70301047

UDI:

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

UBM Calculation

Calculation

Initial date of Tube statistic data log	Jun 25, 2024				
Final date of Tube statistic data log	Jun 25, 2024				
Number of hours					
Total initial count	0				
Total final count	0				
Exposure difference	0				
Number of days	0				
Exposure per year	0				
Categories	High				

PM Start

PM Start

Performing a Visual Inspection	PASS
Switching On the system	PASS
Checking system log files (Checking M Cabinet log files can be skipped if generator errors are not found in system log files)	PASS
Checking system functions	PASS
Checking the thermal overload switch	PASS
Safety, Display and Indicating systems	PASS

	Min	Max	UoM	Measurement	Result	Comments
PM Start						
PM Start						
System Maintenance record					PASS	
Functional Status					PASS	
PM Table TH						
Table TH						
Functional Status					PASS	
Table Base					PASS	
Bucky Unit					PASS	
Final Steps					PASS	
PM Wall Stand VS						
Wall Stand VS (Question set 1)						
Fixations					PASS	
Wall Stand VS (Question set 2)						
Functional Status					PASS	
Functional test vertical movement					PASS	
Covers					PASS	
Vertical Chain					PASS	
Guide Rails and Castors					PASS	
Vertical Collision Check					PASS	
Bucky Unit					PASS	
Options					PASS	
Wrap Up					PASS	
PM Ceiling Suspension						
Functional Status						
Pre checks					PASS	
Functional Status1						
a. Functional Test Long Movement					PASS	
b. Functional Test Lateral Movement					PASS	
c. Functional Test Vertical Movement					PASS	
d. Functional Test Tilt and Swing					PASS	
Test of Geometry						
Vertical test					PASS	
Horizontal test(if wall stand is installed)					PASS	

	Min	Max	UoM	Measurement	Result	Comments
PM Ceiling Suspension						
Covers						
Rear Cover of telescopic carriage					PASS	
Centre cover of telescopic carriage					PASS	
Tube cover					PASS	
Fixations						
Fixations					PASS	
CSM/CSM3 - Component Integrity						
Component integrity					PASS	
Cleaning brush					Pass	
Catches					Pass	
Alpha Drive (CSM3)					Pass	
Cable Roller					Pass	
Cable Carriages					Pass	
X ray tube and display installation					Pass	
CSM/CSM3 - Castors						
Castors					PASS	
CSM/CSM3						
Vertical Column					PASS	
Tracking Unit					PASS	
CSM/CSM3 - Wire Ropes						
Wire Ropes					PASS	
Wrap Up						
Wrap up					PASS	
PM Workstation						
AWS - Acquisition Work Spot						
Functional Status					PASS	
Philips Remote Service Functionality Checked					Customer denied PRS	
Cleaning AWS and Accessories					PASS	
Hard Disk Maintenance					PASS	
Service Partition Maintenance					PASS	
Patient Database Maintenance					PASS	
Updating Virus Guard (if applicable)					PASS	
Updating the Eleva OS SPxx (if applicable)					PASS	

	Min	Max	UoM	Measurement	Result	Comments
PM M-Cabinet						
M-Cabinet (Question set 1)						
Functional Status					PASS	
M-Cabinet (Question set 2)						
Switching off the system (skip if GEST errors like arcing, etc are not present)					PASS	
Turning off the mains power (skip if GEST errors like arcing, etc are not present)					PASS	
Removing the covers (skip if GEST errors like arcing, etc are not present)					PASS	
Cable connections (skip if GEST errors like arcing, etc are not present)					PASS	
M-Cabinet (Question set 3)						
Functional Status					PASS	
Tube calibration					PASS	
Backup of Generator data					PASS	
Final steps					PASS	
PM X-ray Tubes						
SRO/RO/SRM X-ray Tube Assemblies						
Preparatory Work					PASS	
Mechanical Status X-Ray Tube Assy.					PASS	
Status of Cables and Connectors					PASS	
Final Work					PASS	
PM Collimators NICOL/Galileo/XD						
Collimators NICOL/Galileo/XD (Question set 1)						
Collimator Lamp Replacement					PASS	
Safety and Radiation Safety Check					PASS	
Functional Check					PASS	
Collimators NICOL/Galileo/XD (Question set 2)						
General Inspection					PASS	
Cleaning and Lubrication					PASS	
PM Accessories						
Accessories						
Additional Grids					PASS	

	Min	Max	UoM	Measurement	Result	Comments
PM Accessories						
Accessories						
Stitching Support					PASS	
ES Table TH						
Protective Earth Table TH						
2) Potential Balance Point SX80		300	mOhm	92	Pass	
4) Foot Switch		300	mOhm	62	Pass	
5) Transversal carriage at SAYX 13		300	mOhm	100	Pass	
7) Service Cover		300	mOhm	68	Pass	
Protective Earth Bucky unit Table (Fixed Detector)						
2) Grid ejector		300	mOhm	62	Pass	
5) Rear cover screw		300	mOhm	70	Pass	
Touch Current Table TH						
2) Potential Balance Point SX80		20	μA	2.8	Pass	
3) Table Top Side Rail (Front)		20	μA	2.7	Pass	
3) Table Top Side Rail (Back)		20	μA	2.7	Pass	
4) Foot Switch		20	μA	2.7	Pass	
5) Transversal carriage at SAYX 13		20	μA	2.7	Pass	
7) Service Cover		20	μA	2.6	Pass	
Touch Current Bucky unit Table (Fixed Detector)						
2) Grid ejector		20	μA	2.3	Pass	
5) Rear cover screw		20	μA	2.3	Pass	
ES Wall Stand VS						
Protective Earth Wall Stand VS						
2) Column screw		300	mOhm	66	Pass	
3) Wall box cover screw		300	mOhm	90	Pass	
Touch Current Wallstand VS						
2) Column screw		20	μA	2.7	Pass	
3) Wall box cover screw		20	μA	2.7	Pass	
4) Rear cover screw wall stand display		20	μA	2.7	Pass	
5) Holder screw remote control station		20	μA	2.7	Pass	

	Min	Max	UoM	Measurement	Result	Comments
ES Bucky Units Wall Stands						
Protective Earth Bucky unit Wall Stand (Fixed Detector)						
2) Grid ejector		300	mOhm	90	Pass	
6) Rear cover screw		300	mOhm	80	Pass	
Touch Current Bucky unit Wall Stand (Fixed Detector)						
2) Grid ejector		20	µA	2.3	Pass	
6) Rear cover screw		20	µA	2.3	Pass	
ES Ceiling Suspension						
Protective Earth CS						
4) Union nut of the high voltage cable		300	mOhm	60	Pass	
6) Accessory rail screw		300	mOhm	266	Pass	
7) Rear screw of the tube control handle		300	mOhm	76	Pass	
Touch Current CS						
4) Union nut of the high voltage cable		20	µA	2.2	Pass	
5) Key switch		20	µA	2.2	Pass	
6) Accessory rail screw		20	µA	2.2	Pass	
7) Rear screw of the tube control handle		20	µA	2.2	Pass	
8) Clutch disk		20	µA	2.2	Pass	
ES M-Cabinet						
Protective Earth M-Cabinet CXA						
2/B) Front cover		300	mOhm	21	Pass	
3/C) Top cover		300	mOhm	30	Pass	
4/D) Rear cover		300	mOhm	22	Pass	
5/E) Left cover		300	mOhm	30	Pass	
6/F) Right cover		300	mOhm	26	Pass	
Touch Current M-Cabinet CXA						
2/B) Front cover		20	µA	1.5	Pass	
3/C) Top cover		20	µA	1.2	Pass	
4/D) Rear cover		20	µA	1.2	Pass	
5/E) Left cover		20	µA	1.4	Pass	
6/F) Right cover		20	µA	1.5	Pass	
ES Workstation						
Protective Earth Acquisition Workspot						
2) Screw (rear side of PC)		300	mOhm	80	Pass	

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

ES Wall Connection Box

Protective Earth Wall Connection Box

2) Upper Left Bezel	300	mOhm	30	Pass	
3) Lower Left Bezel	300	mOhm	32	Pass	
4) Upper Right Bezel	300	mOhm	36	Pass	
5) Lower Right Bezel	300	mOhm	28	Pass	
6) Front Cover Screw	300	mOhm	30	Pass	

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Vario: 7

Large: 47,528

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Philips Representative Signature:

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Engineer Name: TOMOKI ABE

Date: Jun 25, 2024